

Anti-Collagen VI/Col6A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P12109

Item Number: CM0005481



Introduction

Recombinant anti-Collagen VI/COL6A1 rabbit monoclonal antibody validated for western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat COL6A1 (P12109). Ideal for extracellular matrix and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Collagen VI/COL6A1 Rabbit mAb
Remarks/Alias	OPLL; BTHLM1; UCHMD1; Collagen VI/COL6A1
Species	Human
Gene ID (Human)	1291
Gene ID	1291
Immunogen	Recombinant fusion protein containing amino acids 37-235 of human Collagen VI/COL6A1 (P12109)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12109
Protein Weight	109kDa
Shipping	Ice bag