


Anti-Collagen X/Col10A1 Rabbit Monoclonal Antibody For Wb, Elisa - Q03692

Item Number: CM0002802



Introduction

High-quality rabbit monoclonal antibody targeting human Collagen X/COL10A1 for WB and ELISA applications. Validated in human, mouse, and rat samples. Ideal for cartilage mineralization and skeletal disorder research.

[Learn More](#)

Parameter	Value
Product Name	Collagen X/COL10A1 Rabbit mAb
Remarks/Alias	COL10A1; collagen alpha-1(X) chain; Collagen X/COL10A1
Species	Human
Gene ID (Human)	1300
Gene ID	1300
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 500-600 of human Collagen X/COL10A1 (Q03692).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q03692
Protein Weight	66kDa
Shipping	Ice bag