


Anti-Collagen X/Col10A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q03692

Item Number: CM0022275



Introduction

Collagen X/COL10A1 rabbit polyclonal antibody validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Reacts with human, mouse, rat. Recombinant immunogen. Detects a 66 kDa protein. Suitable for cartilage and bone mineralization research.

[Learn More](#)

Parameter	Value
Product Name	Collagen X/COL10A1 Rabbit pAb
Remarks/Alias	COL10A1; collagen alpha-1(X) chain; Collagen X/COL10A1
Species	Human
Gene ID (Human)	1300
Gene ID	1300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 521-680 of human Collagen X/COL10A1 (NP_000484.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03692
Protein Weight	66 kDa
Shipping	Ice bag