

## Anti-Collagen Xvii/Col17A1 Rabbit Mab For Wb, Ihc-P, Elisa - Q9Umd9

Item Number: CM0004736



### Introduction

Rabbit monoclonal antibody targeting human Collagen XVII (COL17A1, BP180). Validated for WB, IHC-P, ELISA. Ideal for hemidesmosome, basement membrane, and bullous pemphigoid research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Collagen XVII/COL17A1 Rabbit mAb                                                                                                        |
| Remarks/Alias    | ERED; JEB4; BP180; BPA-2; BPAG2; LAD-1; BA16H23.2; Collagen XVII/COL17A1                                                                |
| Species          | Human                                                                                                                                   |
| Gene ID (Human)  | 1308                                                                                                                                    |
| Gene ID          | 1308                                                                                                                                    |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1300-1400 of human Collagen XVII/COL17A1 (Q9UMD9). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                                   |
| Application      | WB, IHC-P, ELISA                                                                                                                        |
| Cross Reactivity | Human                                                                                                                                   |
| SWISS            | Q9UMD9                                                                                                                                  |
| Protein Weight   | 150kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                 |