


Anti-Commd1 Polyclonal Antibody For If/Icc, Elisa - Q8N668

Item Number: CM0027054



Introduction

Rabbit polyclonal antibody against COMMD1 (MURR1) for IF/ICC and ELISA. Detects human, mouse, and rat COMMD1 (21kDa). Suitable for endosomal recycling, NF-kB, and copper homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	COMMD1 Rabbit pAb
Remarks/Alias	MURR1; C2orf5; COMMD1
Species	Human
GeneID (Human)	150684
GeneID	150684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human COMMD1 (NP_689729.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N668
Protein Weight	21kDa
Shipping	Ice bag