


Anti-Commd5 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Gzq3

Item Number: CM0020956



Introduction

Rabbit polyclonal antibody against COMMD5 (HCARG), validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets a scaffold protein involved in endosomal recycling and cell proliferation regulation.

[Learn More](#)

Parameter	Value
Product Name	COMMD5 Rabbit pAb
Remarks/Alias	HCARG; HT002; COMMD5
Species	Human
Gene ID	28991
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human COMMD5 (NP_054785.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ3
Protein Weight	25kDa
Shipping	Ice bag