


Anti-Comp Polyclonal Antibody For Wb, If-P, Elisa - P49747

Item Number: CM0029287



Introduction

Rabbit polyclonal antibody targeting human COMP (Cartilage Oligomeric Matrix Protein), validated for WB, IF-P, ELISA. Cross-reacts with mouse, rat. Used in cartilage ECM and osteoarthritis research.

[Learn More](#)

Parameter	Value
Product Name	COMP Rabbit pAb
Remarks/Alias	MED; CTS2; EDM1; EPD1; TSP5; PSACH; THBS5; TSP-5; COMP
Species	Human
GeneID (Human)	1311
GeneID	1311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-160 of human COMP (NP_000086.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49747
Protein Weight	83 kDa
Shipping	Ice bag