


Anti-Comp Polyclonal Antibody For Wb, If/Icc, Elisa - P49747

Item Number: CM0017455



Introduction

Rabbit polyclonal antibody against human COMP (Thrombospondin-5) for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for cartilage, ECM, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	COMP Rabbit pAb
Remarks/Alias	MED; CTS2; EDM1; EPD1; TSP5; PSACH; THBS5; TSP-5; COMP
Species	Human
GeneID	1311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-160 of human COMP (NP_000086.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49747
Protein Weight	83kDa
Shipping	Ice bag