


Anti-Comp Rabbit Monoclonal Antibody For Wb, Elisa - P49747

Item Number: CM0003526



Introduction

Rabbit monoclonal antibody against human Cartilage Oligomeric Matrix Protein (COMP, TSP5). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for cartilage and ECM research.

[Learn More](#)

Parameter	Value
Product Name	COMP Rabbit mAb
Remarks/Alias	MED; CTS2; EDM1; EPD1; TSP5; PSACH; THB55; TSP-5; COMP
Species	Human
GeneID (human)	1311
GeneID	1311
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human COMP (NP_000086.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P49747
Protein Weight	83 kDa
Shipping	Ice bag