

Anti-Complement C2 Rabbit Polyclonal Antibody For Wb And Elisa - P06681

Item Number: CM0017903



Introduction

Rabbit polyclonal antibody targeting human Complement C2, suitable for Western blot and ELISA. Enables study of complement activation, immune defense, and age-related macular degeneration. Cross-reacts with human C2 (83 kDa).

[Learn More](#)

Parameter	Value
Product Name	C2 Rabbit pAb
Remarks/Alias	CO2; ARMD14; C2
Species	Human
GeneID (Human)	717
GeneID	717
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 453-752 of human C2 (NP_000054.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06681
Protein Weight	83kDa
Shipping	Ice bag