

Anti-Complement C3 Rabbit Polymab For Wb And Elisa - P01024 / P01027

Item Number: CM0009045



Introduction

Rabbit monoclonal antibody against Complement C3 (C3, ASP, C3a, C3b) validated for WB and ELISA. Reacts with human, mouse, and rat C3. Ideal for complement cascade research, inflammation studies, and metabolism investigations.

[Learn More](#)

Parameter	Value
Product Name	Complement C3 Rabbit PolymAb®
Remarks/Alias	ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p
Species	Human/Mouse
GeneID (Human)	718
GeneID	718/12266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1364-1663 of human/mouse Complement C3 (NP_033908.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01024/P01027
Protein Weight	186kDa
Shipping	Ice bag