


Anti-Complement C4A Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P0C0L4

Item Number: CM0010197



Introduction

High-quality rabbit monoclonal antibody against human Complement C4A (SWISS: P0C0L4). Validated for Western blot, immunofluorescence (IF-P), and ELISA applications. Ideal for studying the classical complement pathway and innate immunity.

[Learn More](#)

Parameter	Specification
Product Name	Complement C4A Rabbit mAb
Remarks/Alias	C4; RG; C4S; C04; C4A2; C4A3; C4A4; C4A6; C4AD; CPAMD2; Complement C4A
Species	Human
Gene ID (Human)	720
Gene ID	720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 940-1237 of human Complement C4A (P0C0L4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P0C0L4
Protein Weight	193kDa
Shipping	Ice bag