

Anti-Complement Factor B (Cfb) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P00751

Item Number: CM0007171



Introduction

Rabbit monoclonal antibody against human Complement factor B (CFB). Validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for complement alternative pathway studies. Molecular weight ~86kDa.

[Learn More](#)

Parameter	Value
Product Name	Complement factor B (CFB) Rabbit mAb
Remarks/Alias	BF; FB; BFD; GBG; CFAB; CFBD; PBF2; AHUS4; FBI12; H2-Bf; ARMD14; Complement factor B (CFB)
Species	Human
GeneID (Human)	629
GeneID	629
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human Complement factor B (CFB) (NP_001701.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00751
Protein Weight	86kDa
Shipping	Ice bag