

Anti-Complement Factor B (Cfb) Rabbit Polyclonal Antibody For Wb, Elisa - P00751

Item Number: CM0017240



Introduction

High-quality rabbit polyclonal antibody targeting human Complement Factor B (CFB), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for complement pathway research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Complement factor B (CFB) Rabbit pAb
Remarks/Alias	BF; FB; BFD; GBG; CFAB; CFBD; PBF2; AHUS4; FBI12; H2-Bf; ARMD14; Complement factor B (CFB)
Species	Human
GeneID (Human)	629
GeneID	629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-420 of human Complement factor B (CFB) (NP_001701.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P00751
Protein Weight	86kDa
Shipping	Ice bag