

Anti-Complement Factor B (Cfb) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00751

Item Number: CM0020126



Introduction

Rabbit polyclonal antibody against human Complement factor B (CFB), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying alternative complement pathway activation.

[Learn More](#)

Parameter	Value
Product Name	Complement factor B (CFB) Rabbit pAb
Remarks/Alias	BF; FB; BFD; GBG; CFAB; CFBD; PBF2; AHUS4; FBI12; H2-Bf; ARMD14; Complement factor B (CFB)
Species	Human
Gene ID (Human)	629
Gene ID	629
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00751
Protein Weight	86kDa
Shipping	Ice bag