

Anti-Complement Factor B Polyclonal Antibody For Wb, If/Icc, Elisa - P00751

Item Number: CM0021491



Introduction

Rabbit polyclonal antibody against human Complement factor B (CFB, P00751). Suitable for WB, IF/ICC, ELISA; cross-reacts with human and mouse. A key reagent for alternative complement pathway, innate immunity, and related disease research.

[Learn More](#)

Parameter	Value
Product Name	Complement factor B Rabbit pAb
Remarks/Alias	BF; FB; BFD; GBG; CFAB; CFBD; PBF2; AHUS4; FBI12; H2-Bf; ARMD14; Complement factor B (CFB)
Species	Human
GeneID (human)	629
GeneID	629
Immunogen	Recombinant protein corresponding to amino acids 26-764 of human Complement factor B (NP_001701.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00751
Protein Weight	86kDa
Shipping	Ice bag