

Anti-Complement Factor D Rabbit Monoclonal Antibody For Wb, Elisa - P00746

Item Number: CM0007445



Introduction

High-quality rabbit anti-human Complement factor D (Adipsin) monoclonal antibody validated for Western blot and ELISA applications. Targets the secreted serine protease with a protein weight of 27kDa. Ideal for alternative complement pathway research.

[Learn More](#)

Parameter	Value
Product Name	Complement factor D Rabbit mAb
Remarks/Alias	DF; ADN; PFD; ADIPSIN; Complement factor D
Species	Human
GeneID (Human)	1675
GeneID	1675
Immunogen	Recombinant protein of human Complement factor D
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00746
Protein Weight	27kDa
Shipping	Ice bag