

Anti-Complement Factor D Rabbit Polyclonal Antibody For Wb, Elisa - P00746

Item Number: CM0023486



Introduction

Rabbit polyclonal antibody targeting human Complement factor D (Adipsin/CFD). Validated for WB and ELISA, suitable for studying the alternative complement pathway, innate immunity, and inflammatory diseases. 27 kDa secreted serine protease.

[Learn More](#)

Parameter	Value
Product Name	Complement factor D Rabbit pAb
Remarks/Alias	DF; ADN; PFD; ADIPSIN; Complement factor D
Species	Human
GeneID (Human)	1675
GeneID	1675
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-253 of human Complement factor D (NP_001919.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00746
Protein Weight	27kDa
Shipping	Ice bag