

Anti-Complement Factor H Rabbit Monoclonal Antibody For Wb, Elisa - P08603

Item Number: CM0009954



Introduction

Rabbit monoclonal antibody targeting human complement factor H (CFH), validated for WB and ELISA. Molecular weight 139kDa. Ideal for studying complement regulation, innate immunity, and related diseases.

[Learn More](#)

Parameter	Value
Product Name	Complement factor H Rabbit mAb
Remarks/Alias	FH; HF; HF1; HF2; HUS; FHL1; AHUS1; AMBP1; ARMD4; ARMS1; CFHL3; Complement factor H
Species	Human
Gene ID (Human)	3075
Gene ID	3075
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Complement factor H (P08603).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08603
Protein Weight	139kDa
Shipping	Ice bag