

Anti-Complement Factor H Rabbit Pab For Wb - P08603

Item Number: CM0017931



Introduction

Rabbit polyclonal antibody targeting Complement factor H, a key soluble inhibitor of the complement alternative pathway. Validated for WB, IF/ICC, IF-P, and ELISA in human, mouse, and rat. Predicted molecular weight 139kDa.

[Learn More](#)

Parameter	Value
Product Name	Complement factor H Rabbit pAb
Remarks/Alias	FH; HF; HF1; HF2; HUS; FHL1; AHUS1; AMBP1; ARMD4; ARMS1; CFHL3; Complement factor H
Species	Human
Gene ID (Human)	3075
Gene ID	3075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human Complement factor H (NP_004485.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08603
Protein Weight	139kDa
Shipping	Ice bag