

Anti-Complement Factor H Rabbit Polyclonal Antibody For Wb, IHC-P, IF/ICC, Elisa - P08603

Item Number: CM0022488



Introduction

Rabbit polyclonal antibody against human Complement factor H (CFH), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying complement alternative pathway regulation and age-related macular degeneration.

[Learn More](#)

Parameter	Value
Product Name	Complement factor H Rabbit pAb
Remarks/Alias	FH; HF; HF1; HF2; HUS; FHL1; AHUS1; AMBP1; ARMD4; ARMS1; CFHL3; Complement factor H
Species	Human
Gene ID (Human)	3075
Gene ID	3075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-323 of human Complement factor H (NP_000177.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08603
Protein Weight	139kDa
Shipping	Ice bag