


Anti-Comt Rabbit Monoclonal Antibody For Wb And Elisa - P21964

Item Number: CM0010150



Introduction

Rabbit monoclonal antibody targeting human COMT (catechol O-methyltransferase), suitable for Western blot and ELISA. Cross-reacts with human and rat. Predicted molecular weight 30 kDa. Ideal for neurotransmitter metabolism research.

[Learn More](#)

Parameter	Value
Product Name	COMT Rabbit mAb
Remarks/Alias	HEL-S-98n; COMT
Species	Human
Gene ID (Human)	1312
Gene ID	1312
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human COMT (P21964).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P21964
Protein Weight	30kDa
Shipping	Ice bag