


Anti-Comt Rabbit Polyclonal Antibody For Wb, Elisa - P21964

Item Number: CM0027828



Introduction

Rabbit polyclonal antibody targeting human COMT (Catechol O-methyltransferase), reactive with human, mouse, and rat. Validated for WB and ELISA. Ideal for catecholamine metabolism and neurotransmitter degradation research.

[Learn More](#)

Parameter	Value
Product Name	COMT Rabbit pAb
Remarks/Alias	HEL-S-98n; COMT
Species	Human
GeneID (Human)	1312
GeneID	1312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-221 of human COMT (NP_009294.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21964
Protein Weight	30kDa
Shipping	Ice bag