


Anti-Comt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P21964

Item Number: CM0021102



Introduction

Rabbit polyclonal antibody against catechol O-methyltransferase (COMT) validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for neurotransmitter metabolism studies - P21964.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] COMT Rabbit pAb
Remarks/Alias	HEL-S-98n; [KD Validated] COMT
Species	Human
GeneID (Human)	1312
GeneID	1312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-271 of human COMT (NP_000745.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21964
Protein Weight	30kDa
Shipping	Ice bag