


Anti-Connexin 43 Monoclonal Antibody For Wb And Elisa - P17302

Item Number: CM0007044



Introduction

Rabbit monoclonal Connexin 43 antibody validated for WB, ELISA. Reacts with human, mouse, rat. Targets gap junction alpha-1 protein, key for intercellular signaling and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Connexin 43 Rabbit mAb
Remarks/Alias	HSS; CMDR; CX43; EKVP; GJAL; ODDD; AVSD3; EKVP3; HLHS1; PPKCA; Connexin 43
Species	Human
GeneID (Human)	2697
GeneID	2697
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-382 of human Connexin 43 (NP_000156.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17302
Protein Weight	43 kDa
Shipping	Ice bag