


Anti-Connexin 43 Polyclonal Antibody For Wb, Elisa - P17302

Item Number: CM0023271



Introduction

Rabbit polyclonal antibody targeting human Connexin 43 (Cx43/GJA1). Suitable for WB and ELISA, cross-reactive with mouse and rat. Ideal for gap junction biology and antiviral innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Connexin 43 Rabbit pAb
Remarks/Alias	HSS; CMDR; CX43; EKVP; GJAL; ODDD; AVSD3; EKVP3; HLHS1; PPKCA; Connexin 43
Species	Human
GeneID (Human)	2697
GeneID	2697
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 283-382 of human Connexin 43 (NP_000156.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17302
Protein Weight	43kDa
Shipping	Ice bag