


## Anti-Connexin 43 Rabbit Polymab® For Wb, If/Icc, If-P, Elisa - P17302

Item Number: CM0006324



### Introduction

Connexin 43 Rabbit PolymAb® (P17302) for human, mouse, rat WB, IF/ICC, IF-P, ELISA. High-quality monoclonal antibody for gap junction and innate immunity research. Reliable IVD raw material.

[Learn More](#)

| Parameter        | Value                                                                                                                            |
|------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Connexin 43 Rabbit PolymAb®                                                                                                      |
| Remarks/Alias    | HSS; CMDR; CX43; EKVP; GJAL; ODDD; AVSD3; EKVP3; HLHS1; PPKCA                                                                    |
| Species          | Human                                                                                                                            |
| GeneID (human)   | 2697                                                                                                                             |
| GeneID           | 2697                                                                                                                             |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-382 of human Connexin 43 (NP_000156.1). |
| Source           | Rabbit                                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                                            |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                |
| SWISS            | P17302                                                                                                                           |
| Protein Weight   | 43kDa                                                                                                                            |
| Shipping         | Ice bag                                                                                                                          |