

## Anti-Contactin 2 Monoclonal Antibody For Wb And Elisa - Q02246

Item Number: CM0010416



### Introduction

Rabbit monoclonal antibody against human Contactin-2, validated for Western blot and ELISA, cross-reacts with mouse. Suitable for neuroscience research on axonal domains and cell adhesion.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Contactin 2 Rabbit mAb                                                                                                      |
| Remarks/Alias    | AXT; TAX; TAX1; FAME5; TAG-1; Contactin 2                                                                                   |
| Species          | Human                                                                                                                       |
| Gene ID (Human)  | 6900                                                                                                                        |
| Gene ID          | 6900                                                                                                                        |
| Immunogen        | Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Contactin 2 (Q02246). |
| Source           | Rabbit                                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                                       |
| Application      | WB, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                                |
| SWISS            | Q02246                                                                                                                      |
| Protein Weight   | 113kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                     |