


Anti-Contactin 2 Rabbit Polyclonal Antibody For Wb, Elisa - Q02246

Item Number: CM0027589



Introduction

Rabbit polyclonal antibody against human Contactin 2 (TAG-1/TAX-1). Suitable for WB and ELISA, cross-reacts with mouse and rat. Target involved in axonal domain organization and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	Contactin 2 Rabbit pAb
Remarks/Alias	AXT; TAX; TAX1; FAME5; TAG-1; Contactin 2
Species	Human
Gene ID (human)	6900
Gene ID	6900
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 785-1012 of human Contactin 2 (NP_005067.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q02246
Protein Weight	113kDa
Shipping	Ice bag