

Anti-Copa Monoclonal Antibody For Wb, Ihc-P, Elisa - P53621

Item Number: CM0003629



Introduction

Recombinant monoclonal antibody targeting COPA (alpha-COP, coatmer subunit alpha). Validated for Western blot, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for intracellular trafficking and COPI vesicle research.

[Learn More](#)

Parameter	Value
Product Name	COPA Rabbit mAb
Remarks/Alias	AILJK; HEP-COP; alpha-COP; COPA
Species	Human
GeneID (Human)	1314
GeneID	1314
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human COPA (P53621).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53621
Protein Weight	138kDa
Shipping	Ice bag