

Anti-Copb2 Rabbit Monoclonal Antibody For Wb, Elisa - P35606

Item Number: CM0002348

Introduction

Rabbit monoclonal antibody against COPB2, validated for Western blot and ELISA. Detects human and mouse COPB2 (coatamer subunit beta') involved in ER-Golgi transport, 102kDa. Ideal for research in trafficking and related diseases.

[Learn More](#)

Parameter	Value
Product Name	COPB2 Rabbit mAb
Remarks/Alias	OPDD; MCPH19; beta'-COP; COPB2
Species	Human
Gene ID (Human)	9276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 657-906 of human COPB2 (NP_004757.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35606
Protein Weight	102kDa
Shipping	Ice bag