


Anti-Copb2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P35606

Item Number: CM0006747



Introduction

Anti-COPB2 rabbit monoclonal antibody validated for WB, IF/ICC, ELISA in human samples. Targets coatamer subunit beta'-COP (p102), a key component of COPI vesicles important for Golgi-to-ER retrograde transport.

[Learn More](#)

Parameter	Value
Product Name	COPB2 Rabbit mAb
Remarks/Alias	OPDD; MCPH19; beta'-COP; COPB2
Species	Human
GeneID (Human)	9276
GeneID	9276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 657-906 of human COPB2 (NP_004757.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P35606
Protein Weight	102kDa
Shipping	Ice bag