


Anti-Copb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35606

Item Number: CM0030066



Introduction

Rabbit polyclonal antibody targeting COPB2, suitable for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat COPB2 (~102 kDa). Ideal for studying coatamer complex function in Golgi trafficking and protein transport.

[Learn More](#)

Parameter	Value
Product Name	COPB2 Rabbit pAb
Remarks/Alias	OPDD; MCPH19; beta'-COP; COPB2
Species	Human
Gene ID (Human)	9276
Gene ID	9276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 657-906 of human COPB2 (NP_004757.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35606
Protein Weight	102kDa
Shipping	Ice bag