


Anti-Cope Rabbit Polyclonal Antibody For Wb And Elisa - O14579

Item Number: CM0017304



Introduction

Rabbit polyclonal antibody targeting human COPE (epsilon-COP), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Useful for studying COPI vesicle transport and Golgi function.

[Learn More](#)

Parameter	Value
Product Name	COPE Rabbit pAb
Remarks/Alias	epsilon-COP; COPE
Species	Human
GeneID (human)	11316
GeneID	11316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-308 of human COPE (NP_009194.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14579
Protein Weight	34kDa
Shipping	Ice bag