


Anti-Copg1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y678

Item Number: CM0017689



Introduction

COPG1 rabbit polyclonal antibody for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, rat COPG1 (coatamer subunit gamma-1, ~98 kDa). Reliable for Golgi transport and COPI vesicle research.

[Learn More](#)

Parameter	Value
Product Name	COPG1 Rabbit pAb
Remarks/Alias	COPG; COPG1
Species	Human
Gene ID (Human)	22820
Gene ID	22820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 555-874 of human COPG1 (NP_057212.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y678
Protein Weight	98kDa
Shipping	Ice bag