

Anti-Copg2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Ubf2

Item Number: CM0016527

Introduction

High-quality rabbit polyclonal antibody targeting COPG2 (gamma-2-COP), a coatamer subunit involved in Golgi vesicle transport. Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat samples. Ideal for studies on retrograde transport and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	COPG2 Rabbit pAb
Remarks/Alias	2-COP; gamma-2-COP; COPG2
Species	Human
GeneID (Human)	26958
GeneID	26958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 542-871 of human COPG2 (NP_036265.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBF2
Protein Weight	98kDa
Shipping	Ice bag