

Anti-Cops2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61201

Item Number: CM0027662



Introduction

Rabbit polyclonal antibody against human COPS2 (COP9 signalosome subunit 2) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat COPS2, a key regulator of ubiquitin-mediated proteolysis. Ideal for cancer and neuronal differentiation research.

[Learn More](#)

Parameter	Value
Product Name	COPS2 Rabbit pAb
Remarks/Alias	CSN2; SGN2; ALIEN; TRIP15; COPS2
Species	Human
Gene ID (Human)	9318
Gene ID	9318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human COPS2 (NP_004227.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61201
Protein Weight	52kDa
Shipping	Ice bag