


Anti-Cops3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uns2

Item Number: CM0030559



Introduction

Rabbit polyclonal antibody targeting human COPS3, a COP9 signalosome subunit regulating ubiquitin-dependent degradation. Suitable for WB, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for signalosome and cullin-RING ligase studies.

[Learn More](#)

Parameter	Value
Product Name	COPS3 Rabbit pAb
Remarks/Alias	CSN3; SGN3; COPS3
Species	Human
GeneID (Human)	8533
GeneID	8533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-423 of human COPS3 (NP_003644.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNS2
Protein Weight	48kDa
Shipping	Ice bag