

Anti-Cops3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Uns2

Item Number: CM0002684



Introduction

Recombinant rabbit monoclonal antibody against human COPS3 (COP9 signalosome complex subunit 3). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for ubiquitin pathway and signalosome research.

[Learn More](#)

Parameter	Value
Product Name	COPS3 Rabbit mAb
Remarks/Alias	CSN3; SGN3; COPS3
Species	Human
GeneID (Human)	8533
GeneID	8533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-423 of human COPS3 (Q9UNS2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNS2
Protein Weight	48kDa
Shipping	Ice bag