


Anti-Cops6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7L5N1

Item Number: CM0029756



Introduction

COPS6 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA against human, mouse, rat; 36kDa COP9 signalosome subunit 6. Ideal for ubiquitin pathway, signalosome, and deneddylation research.

[Learn More](#)

Specification	Value
Product Name	COPS6 Rabbit pAb
Remarks/Alias	CSN6; MOV34-34KD; COPS6
Species	Human
GeneID (Human)	10980
GeneID	10980
Immunogen	Recombinant fusion protein containing amino acids 58-327 of human COPS6 (NP_006824.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L5N1
Protein Weight	36kDa
Shipping	Ice bag