


Anti-Cops8 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99627

Item Number: CM0015522



Introduction

COPS8 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects COPS8 in human, mouse, and rat. Targets COP9 signalosome subunit 8 for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	COPS8 Rabbit pAb
Remarks/Alias	COP9; CSN8; SGN8; COPS8
Species	Human
GeneID (Human)	10920
GeneID	10920
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-209 of human COPS8 (NP_006701.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99627
Protein Weight	23kDa
Shipping	Ice bag