

Anti-Copz1 Polyclonal Antibody For Wb, If, Elisa - P61923

Item Number: CM0013431



Introduction

Coatamer subunit zeta-1 (COPZ1) rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Suitable for human and mouse samples; targets a ~20kDa protein involved in Golgi-to-ER transport. Ideal for molecular biology and IVD research.

[Learn More](#)

Parameter	Specifications
Product Name	COPZ1 Rabbit pAb
Remarks/Alias	COPZ; CGI-120; HSPC181; zeta-COP; zeta1-COP; COPZ1
Species	Human
GeneID (Human)	22818
GeneID	22818
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61923
Protein Weight	20kDa
Shipping	Ice bag