

Anti-Coq10B Rabbit Polyclonal Antibody For Wb And Elisa - Q9H8M1

Item Number: CM0014635

Introduction

Rabbit polyclonal antibody against human COQ10B (Coenzyme Q-binding protein COQ10 homolog B, mitochondrial). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for respiratory chain studies.

[Learn More](#)

Parameter	Value
Product Name	COQ10B Rabbit pAb
Remarks/Alias	COQ10B
Species	Human
GeneID	80219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-150 of human COQ10B (NP_079423.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H8M1
Protein Weight	27kDa
Shipping	Ice bag