


Anti-Coq7 Rabbit Polyclonal Antibody For Wb And Elisa - Q99807

Item Number: CM0027974



Introduction

High-quality rabbit polyclonal antibody targeting human COQ7 (NADPH-dependent 3-demethoxyubiquinone 3-hydroxylase), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for ubiquinone biosynthesis and mitochondrial stress response studies.

[Learn More](#)

Parameter	Value
Product Name	COQ7 Rabbit pAb
Remarks/Alias	CAT5; CLK1; CLK-1; COQ10D8; COQ7
Species	Human
Gene ID	10229
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-217 of human COQ7 (NP_057222.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99807
Protein Weight	24kDa
Shipping	Ice bag