

Anti-Cor28 Polyclonal Antibody For Wb And Elisa - Q8L983

Item Number: CM0021097

Introduction

Rabbit polyclonal antibody against Arabidopsis thaliana COR28 for Western blot and ELISA. Suitable for circadian clock and cold stress studies in plant research.

[Learn More](#)

Parameter	Value
Product Name	COR28 Rabbit pAb
Remarks/Alias	F17I5.170; F17I5_170; COR28
Species	Arabidopsis thaliana
GeneID	829544
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-170 of Arabidopsis thaliana COR28 (NP_567946.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q8L983
Protein Weight	25kDa
Shipping	Ice bag