

Anti-Coro1B Rabbit Polyclonal Antibody For Wb And Elisa - Q9Br76

Item Number: CM0016920



Introduction

Rabbit polyclonal antibody against human CORO1B (Coronin-1B), validated for Western blot and ELISA. Cross-reactive with rat. Ideal for cell motility, actin dynamics, and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	CORO1B Rabbit pAb
Remarks/Alias	CORONIN-2; CORO1B
Species	Human
Gene ID	57175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-489 of human CORO1B (NP_065174.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9BR76
Protein Weight	54kDa
Shipping	Ice bag