

Anti-Cortactin Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q14247

Item Number: CM0008820



Introduction

Recombinant rabbit monoclonal antibody against human Cortactin (CTTN/EMS1), validated for WB, IP, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying actin cytoskeleton and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Cortactin Rabbit mAb
Remarks/Alias	EMS1; Cortactin
Species	Human
GeneID (Human)	2017
GeneID	2017
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cortactin (Q14247)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14247
Protein Weight	62kDa
Shipping	Ice bag