


Anti-Cortactin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14247

Item Number: CM0017992



Introduction

Cortactin Rabbit Polyclonal Antibody detects endogenous cortactin (aa 1-200) in human, mouse, and rat samples via WB, IHC-P, and ELISA. Ideal for investigating actin dynamics, cell migration, and tumor invasiveness.

[Learn More](#)

Parameter	Value
Product Name	Cortactin Rabbit pAb
Remarks/Alias	EMS1; Cortactin
Species	Human
Gene ID (Human)	2017
Gene ID	2017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Cortactin (NP_612632.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14247
Protein Weight	62kDa
Shipping	Ice bag