


Anti-Cotl1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14019

Item Number: CM0029944



Introduction

Rabbit polyclonal antibody against coactosin-like protein (COTL1/CLP), validated for WB, IHC-P, and ELISA in human, mouse, rat. Targets the 16 kDa actin-binding and 5LO chaperone protein. Ideal for cytoskeleton and leukotriene biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	COTL1 Rabbit pAb
Remarks/Alias	CLP; COTL1
Species	Human
Gene ID	23406
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human COTL1 (NP_066972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14019
Protein Weight	16kDa
Shipping	Ice bag