


Anti-Coup-Tfii/Nr2F2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P24468

Item Number: CM0007431



Introduction

Rabbit monoclonal antibody against human COUP-TFII/NR2F2 (ARP-1), validated for WB, IHC-P, IF/ICC, IP, ELISA, cross-reacts with mouse and rat. Ideal for nuclear receptor and transcription regulation research.

[Learn More](#)

Specification	Details
Product Name	COUP-TFII/NR2F2 Rabbit mAb
Remarks/Alias	ARP1; ARP-1; CHTD4; NF-E3; SRXX5; SVP40; COUPTF2; COUPTFB; TFCOUP2; COUPTFII
Species	Human
GeneID (Human)	7026
GeneID	7026
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 5-100 of human COUP-TFII/NR2F2 (NP_066285.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24468
Protein Weight	46kDa
Shipping	Ice bag